



Sustainably Designed Medical Devices and Surgical Products

# NGS35W Skin Stapler



## Clinical Performance

- Audible and tactile feedback when firing cycle completes
- Alignment indicator designed for accurate staple placement
- Ergonomically designed for ease of use
  - Accommodates a variety of hand positions and sizes
  - Low force-to-fire
- Easy to see staple remaining indicator
- 35 surgical stainless steel staples
- Quality assurance - each device functionally tested

## Sustainable Design

- Made with 69% plant-based material\*
- 67% reduction in energy used for production\*
- Less plastic produced and discarded\*
- Lightweight design can reduce transportation and disposal cost\*

\* Based on comparison to like product made with plastic.

## Staple Specifications

# of staples	Size	Closed Position	Wire diameter
35	12.0 X 3.3 mm 	6.9 x 3.9 mm 	0.58 mm

Product Code	Description	Units of Measure
NGS35W	35 Wide Single-Use Skin Stapler <sup>†</sup>	6 Per Box
NGS35W-24		24 Per Case

US Patents 9226749, 9572575


Always refer to the Instructions for Use included with the product for complete indications, instructions, warnings, precautions and contraindications.

<sup>†</sup> Available bulk, non-sterile for procedure packs through authorized distributors.

NewGen Surgical Inc • 41 Simms Street, Suite B • San Rafael, CA 94901 • Phone: 855-295-4500 • [www.newgensurgical.com](http://www.newgensurgical.com)

Smart Sustainable Design™

# Smart Sustainable Design™



By rethinking and redesigning disposable single-use surgical products used every day in the operating room, NewGen Surgical reduces or eliminates plastic waste while still delivering clinical performance at a price that brings value to our customers.

## Contact Us

Telephone: 855-295-4500

Fax: 415-532-2290

Email: [info@newgensurgical.com](mailto:info@newgensurgical.com)

Follow:   @NewGenSurgical

## About NewGen Surgical

NewGen Surgical, Inc. develops and manufactures sustainably designed, single-use medical devices and surgical products. NewGen Surgical's mission is to minimize plastic waste in the O.R. by redesigning disposable products. NewGen Surgical offers hospitals and surgery centers a solution to their sustainability programs and meet Environmentally Preferable Purchasing (EPP) objectives. Founded in 2012 in the San Francisco Bay Area, NewGen Surgical is a member of Practice Greenhealth's Greening the OR® Initiative and was commended in the 2015 Circular Economy Awards.



NewGen Surgical is a proud supporter of Practice Greenhealth's Greening the OR® Initiative. Practice Greenhealth is the nation's leading healthcare community dedicated to empowering its members to increase their efficiencies in environmental stewardship.



Think Beyond Plastic™ is a social-impact venture addressing plastic pollution by harnessing the forces of entrepreneurship, innovation, and the power of global markets. In 2012, NewGen Surgical became one of their first members of the Think Beyond Plastic™ Accelerator; the first Accelerator in the world to focus on plastic pollution as an innovation opportunity.



NewGen Surgical was commended in the 2015 Circular Economy Awards by the judges of The Young Global Leaders Award, an initiative of the World Economic Forum, in the Circular Economy Entrepreneurs category for early stage organizations at the vanguard of the circular economy demonstrating innovation and market disruption.

 Printed on recycled paper. Please recycle.